

## 兔抗 ZNF512 多克隆抗体

中文名称：兔抗 ZNF512 多克隆抗体

英文名称：Anti-ZNF512 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻 (-20℃)

抗 原：ZNF512

宿 主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclona

### 技术规格

<b>Background:</b>	This gene encodes a protein containing four putative zinc finger motifs. Zinc finger motifs may bind to proteins or nucleic acids. Zinc finger-containing proteins are involved in a variety of processes, including regulation of transcription. Alternative splicing results in multiple transcript variants for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	ZNF512
<b>Immunogen:</b>	Fusion protein of human ZNF512
<b>Full name:</b>	zinc finger protein 512
<b>SwissProt:</b>	Q96ME7
<b>ELISA Recommended dilution:</b>	5000-10000

<b>IHC positive control:</b>	Human colorectal cancer and human cervical cancer
<b>IHC Recommend dilution:</b>	25-100



