

兔抗 LAMA4 多克隆抗体

中文名称：兔抗 LAMA4 多克隆抗体

英文名称：Anti-LAMA4 rabbit polyclonal antibody

别名：LAMA3; CMD1JJ; LAMA4*-1

抗原：LAMA4

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human

相关类别：一抗

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

IHC Recommend dilution:	25-100
IHC positive control:	Human thyroid cancer
ELISA Recommended dilution:	2000-5000
SwissProt:	Q16363
Synonyms :	LAMA3; CMD1JJ; LAMA4*-1
Full name:	laminin, alpha 4
Immunogen:	Fusion protein of human LAMA4
Name of antibody:	LAMA4
Applications:	ELISA, IHC
Background:	Laminins, a family of extracellular matrix glycoproteins, are t

The major noncollagenous constituent of basement membranes. They have been implicated in a wide variety of biological processes including cell adhesion, differentiation, migration, signaling, neurite outgrowth and metastasis. Laminins are composed of 3 non identical chains: laminin alpha, beta and gamma (formerly A, B1, and B2, respectively) and they form a cruciform structure consisting of 3 short arms, each formed by a different chain, and a long arm composed of all 3 chains.

