

兔抗 TRIM14 多克隆抗体

中文名称: 兔抗 TRIM14 多克隆抗体

英文名称: Anti-TRIM14 rabbit polyclonal antibody

相关类别: 一抗

储 存: 冷冻 (-20℃)

宿 主: Rabbit

抗 原: TRIM14

反应种属: Human

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. The protein localizes to cytoplasmic bodies and its function has not been determined. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	TRIM14
Immunogen:	Synthetic peptide of human TRIM14
Full name:	tripartite motif containing 14
SwissProt:	Q14142

ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	25-100

