

兔抗 SNAP29 多克隆抗体

中文名称：兔抗 SNAP29 多克隆抗体

英文名称： Anti-SNAP29 rabbit polyclonal antibody

别名： CEDNIK; SNAP-29

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SNAP29

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This gene, a member of the SNAP25 gene family, encodes a protein involved in multiple membrane trafficking steps. Two other members of this gene family, SNAP23 and SNAP25, encode proteins that bind a syntaxin protein and mediate synaptic vesicle membrane docking and fusion to the plasma membrane. The protein encoded by this gene binds tightly to multiple syntaxins and is localized to intracellular membrane structures rather than to the plasma membrane. While the protein is mostly membrane-bound, a significant fraction of it is found free in the cytoplasm. Use of mul

	tiple polyadenylation sites has been noted for this gene.
Applications:	ELISA, IHC
Name of antibody:	SNAP29
Immunogen:	Fusion protein of human SNAP29
Full name:	synaptosomal-associated protein, 29kDa
Synonyms :	CEDNIK; SNAP-29
SwissProt:	O95721
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human liver cancer
IHC Recommend dilution:	30-150

