

兔抗 TRIM16 多克隆抗体

中文名称: 兔抗 TRIM16 多克隆抗体

英文名称: Anti-TRIM16 rabbit polyclonal antibody

别 名: EBBP

相关类别: 一抗

储 存: 冷冻(-20℃)

宿 主: Rabbit

抗 原: TRIM16

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:

This gene was identified as an estrogen and anti-estrogen regulated gene in epithelial cells stably expressing estroge n receptor. The protein encoded by this gene contains tw o B box domains and a coiled-coiled region that are char acteristic of the B box zinc finger protein family. The prot eins of this family have been reported to be involved in a variety of biological processes including cell growth, diff erentiation and pathogenesis. Expression of this gene was detected in most tissues. Its function, however, has not ye t been determined.



Applications:	ELISA, IHC
Name of antibody:	TRIM16
Immunogen:	Fusion protein of human TRIM16
Full name:	tripartite motif containing 16
Synonyms:	EBBP
SwissProt:	O95361
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and human esophagus cancer
IHC Recommend dilution:	25-100

