

兔抗 ZNF560 多克隆抗体

中文名称：兔抗 ZNF560 多克隆抗体

英文名称：Anti-ZNF560 rabbit polyclonal antibody

别名：zinc finger protein 560

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ZNF560

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	May be involved in transcriptional regulation.
Applications:	ELISA, IHC
Name of antibody:	ZNF560
Immunogen:	Synthetic peptide of human ZNF560
Full name:	zinc finger protein 560
SwissProt:	Q96MR9
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human ovarian cancer
IHC Recommend dilution:	30-150

