

兔抗 ZNF580 多克隆抗体

中文名称：兔抗 ZNF580 多克隆抗体

英文名称： Anti-ZNF580 rabbit polyclonal antibody

别名： zinc finger protein 580

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： ZNF580

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Involved in the regulation of endothelial cell proliferation and migration. Mediates H ₂ O ₂ -induced leukocyte chemotaxis by elevating interleukin-8 production and may play a role in inflammation. May be involved in transcriptional regulation.
Applications:	ELISA, IHC
Name of antibody:	ZNF580
Immunogen:	Fusion protein of human ZNF580
Full name:	zinc finger protein 580
SwissProt:	Q9UK33
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human breast cancer
IHC Recommend dilution:	30-150

