

## 兔抗 SNAP91 多克隆抗体

中文名称： 兔抗 SNAP91 多克隆抗体

英文名称： Anti-SNAP91 rabbit polyclonal antibody

别名： CALM; AP180

相关类别： 一抗

抗原： SNAP91

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

<b>Background:</b>	Clathrin coat assembly protein AP180 is a protein that in humans is encoded by the SNAP91 gene. Adaptins are components of the adapter complexes which link clathrin to receptors in coated vesicles. Clathrin-associated protein complexes are believed to interact with the cytoplasmic tails of membrane proteins, leading to their selection and concentration. Binding of AP180 to clathrin triskelion induces their assembly into 60-70 nm coats.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SNAP91
<b>Immunogen:</b>	Fusion protein of human SNAP91
<b>Full name:</b>	synaptosomal-associated protein, 91kDa
<b>Synonyms:</b>	CALM; AP180
<b>SwissProt:</b>	O60641

<b>IHC positive control:</b>	Human gastric cancer and human colon cancer
<b>IHC Recommend dilution:</b>	25-100

