

## 兔抗 LAMP1 多克隆抗体

中文名称：兔抗 LAMP1 多克隆抗体

英文名称：Anti-LAMP1 rabbit polyclonal antibody

别名：LAMPA; CD107a; LGP120

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：LAMP1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a member of a family of membrane glycoproteins. This glycoprotein provides selectins with carbohydrate ligands. It may also play a role in tumor cell metastasis.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	LAMP1
<b>Immunogen:</b>	Synthesized peptide derived from internal of human LAMP 1.
<b>Full name:</b>	lysosomal-associated membrane protein 1
<b>Synonyms:</b>	LAMPA; CD107a; LGP120
<b>SwissProt:</b>	P11279

<b>WB Predicted band size:</b>	45 kDa
<b>WB Positive control:</b>	K562, COLO205, Hela, Jurkat, HepG2, Huvec cells lysates
<b>WB Recommended dilution:</b>	500-3000



