

兔抗 LAMP2 多克隆抗体

中文名称:	兔抗 LAMP2 多克隆抗体
英文名称:	Anti-LAMP2 rabbit polyclonal antibody
别名:	LAMPB; CD107b; LAMP-2; LGP110
抗原:	LAMP2
储存:	冷冻 (-20℃)
宿主:	Rabbit
相关类别:	一抗
反应种属:	Human
标记物:	Unconjugate
克隆类型:	rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of a family of membrane glycoproteins. This glycoprotein provides selectins with carbohydrate ligands. It may play a role in tumor cell metastasis. It may also function in the protection, maintenance, and adhesion of the lysosome. Alternative splicing of this gene results in multiple transcript variants encoding distinct proteins
Applications:	ELISA, WB, IHC
Name of antibody:	LAMP2
Immunogen:	Fusion protein of human LAMP2
Full name:	lysosomal-associated membrane protein 2
Synonyms:	LAMPB; CD107b; LAMP-2; LGP110

SwissProt:	P13473
IHC positive control:	Human gastric cancer and human liver cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	45 kDa
WB Positive control:	Human placenta tissue
WB Recommended dilution:	500-2000



