

## 兔抗 LAMP3 多克隆抗体

中文名称：兔抗 LAMP3 多克隆抗体

英文名称：Anti-LAMP3 rabbit polyclonal antibody

别名：LAMP, CD208, DCLAMP, LAMP-3, TSC403, DC LAMP, DC-LAMP

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：LAMP3

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Dendritic cells (DCs) are the most potent antigen-presenting cells. Immature DCs efficiently capture antigens and differentiate into interdigitating dendritic cells (IDCs) in lymphoid tissues that induce primary T-cell responses.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	LAMP3
<b>Immunogen:</b>	Synthetic peptide of human LAMP3
<b>Full name:</b>	lysosomal associated membrane protein 3
<b>Synonyms:</b>	LAMP; CD208; DCLAMP; LAMP-3; TSC403; DC LAMP; DC-LAMP
<b>SwissProt:</b>	Q9UQV4

<b>ELISA Recommended dilution :</b>	5000-10000
<b>IHC positive control:</b>	Human brain and Human colorectal cancer
<b>IHC Recommend dilution:</b>	50-100
<b>WB Predicted band size:</b>	44 kDa
<b>WB Positive control:</b>	Jurkat and Hela cell lysates
<b>WB Recommended dilution:</b>	500-2000



