

兔抗 SNCG 多克隆抗体

中文名称：兔抗 SNCG 多克隆抗体

英文名称：Anti-SNCG rabbit polyclonal antibody

相关类别：一抗

抗 原：SNCG

储 存：冷冻（-20℃）

宿 主：Rabbit

反应种属：Human

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a member of the synuclein family of proteins which are believed to be involved in the pathogenesis of neurodegenerative diseases. Mutations in this gene have also been associated with breast tumor development.
Applications:	ELISA, IHC
Name of antibody:	SNCG
Immunogen:	Fusion protein of human SNCG
Full name:	synuclein, gamma (breast cancer-specific protein 1)
SwissProt:	O76070
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and Human liver cancer
IHC Recommend dilution:	50-200

