

## 兔抗 ZNF672 多克隆抗体

中文名称：兔抗 ZNF672 多克隆抗体

英文名称：Anti-ZNF672 rabbit polyclonal antibody

相关类别：一抗

抗 原：ZNF672

储 存：冷冻（-20℃）

宿 主：Rabbit

反应种属：Human

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	May be involved in transcriptional regulation.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	ZNF672
<b>Immunogen:</b>	Fusion protein of human ZNF672
<b>Full name:</b>	zinc finger protein 672
<b>SwissProt:</b>	Q499Z4
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human thyroid cancer
<b>IHC Recommend dilution:</b>	100-300
<b>WB Predicted band size:</b>	50 kDa
<b>WB Positive control:</b>	SKOV3 cell lysate
<b>WB Recommended dilution:</b>	500-2000

