

兔抗 TRIM38 多克隆抗体

中文名称：兔抗 TRIM38 多克隆抗体

英文名称： Anti-TRIM38 rabbit polyclonal antibody

别名： RNF15; RORET

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： TRIM38

宿主： Rabbit

反应种属： Human Mouse

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes a member of the tripartite motif (TRIM) family. The encoded protein contains a RING-type zinc finger, B box-type zinc finger and SPRY domain. The function of this protein has not been identified. A pseudogene of this gene is located on the long arm of chromosome 4.
Applications:	WB
Name of antibody:	TRIM38
Immunogen:	Synthesized peptide derived from internal of human TRIM38.
Full name:	tripartite motif containing 38

Synonyms :	RNF15; RORET
SwissProt:	O00635
WB Predicted band size:	53 kDa
WB Positive control:	RAW264.7 cells and LOVO cells lysates
WB Recommended dilution:	500-3000

