

## 兔抗 ZNF695 多克隆抗体

中文名称：兔抗 ZNF695 多克隆抗体

英文名称：Anti-ZNF695 rabbit polyclonal antibody

别名：SBZF3

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：ZNF695

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

### 技术规格

<b>Background:</b>	May be involved in transcriptional regulation.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	ZNF695
<b>Immunogen:</b>	Synthesized peptide derived from C-terminal of human ZNF695.
<b>Full name:</b>	zinc finger protein 695
<b>Synonyms :</b>	SBZF3
<b>SwissProt:</b>	Q8IW36
<b>WB Predicted band size:</b>	60 kDa
<b>WB Positive control:</b>	Jurkat and 293 cells lysates
<b>WB Recommended dilution:</b>	500-3000

