

兔抗 TRIM45 多克隆抗体

中文名称: 兔抗 TRIM45 多克隆抗体

英文名称: Anti-TRIM45 rabbit polyclonal antibody

别名: RNF99

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: TRIM45

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a member of the tripartite motif family. The encoded protein may function as a transcriptional repressor of the mitogen-activated protein kinase pathway. Alternatively spliced transcript variants have been described.
Applications:	ELISA, WB, IHC
Name of antibody:	TRIM45
Immunogen:	Synthetic peptide of human TRIM45
Full name:	tripartite motif containing 45
Synonyms :	RNF99

SwissProt:	Q9H8W5
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and Human thyroid cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	64 kDa
WB Positive control:	231, Jurkat, A549, hela and 293T cells
WB Recommended dilution:	500-2000



