

## 兔抗 TRIM49 多克隆抗体

中文名称：兔抗 TRIM49 多克隆抗体

英文名称： Anti-TRIM49 rabbit polyclonal antibody

别名： RNF18; TRIM49A; TRIM49L2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： TRIM49

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene contains a RING zinc finger, a motif known to be involved in protein-protein interactions. This gene has been found to be preferentially expressed in testis. Related pseudogenes and gene duplications have also been identified on chromosome 11.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	TRIM49
<b>Immunogen:</b>	Fusion protein of human TRIM49
<b>Full name:</b>	tripartite motif containing 49
<b>Synonyms :</b>	RNF18; TRIM49A; TRIM49L2

<b>SwissProt:</b>	P0CI25
<b>ELISA Recommended dilution:</b>	2000-10000
<b>IHC positive control:</b>	Human thyroid cancer
<b>IHC Recommend dilution:</b>	30-150

