

兔抗 SNX5 多克隆抗体

中文名称：兔抗 SNX5 多克隆抗体

英文名称： Anti-SNX5 rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻 (-20℃)

宿 主： Rabbit

抗 原： SNX5

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the sorting nexin family. Members of this family contain a phosphoinositide binding domain, which is a phosphoinositide binding domain, and are involved in intracellular trafficking. This protein functions in endosomal sorting, the phosphoinositide signaling pathway, and macropinocytosis. This gene may play a role in the tumorigenesis of papillary thyroid carcinoma. Alternative splicing results in multiple transcript variants encoding different isoforms.
Applications:	ELISA, WB, IHC
Name of antibody:	SNX5
Immunogen:	Synthetic peptide of human SNX5
Full name:	sorting nexin 5

SwissProt:	Q9Y5X3
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human lung cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	47 kDa
WB Positive control:	K562 and Jurkat cells
WB Recommended dilution:	500-2000



