

## 兔抗 SNX9 多克隆抗体

中文名称: 兔抗 SNX9 多克隆抗体

英文名称: Anti-SNX9 rabbit polyclonal antibody

别 名: SNX9; MST155; MSTP155; SDP1; SH3PX1; SH3PXD3A; WISP

相关类别: 一抗

储 存: 冷冻(-20℃) 避光

宿 主: Rabbit

抗 原: SNX9

反应种属: Human, Mouse, Rat

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

## 技术规格

**Background:** 

This gene encodes a member of the sorting ne xin family. Members of this family contain a ph ox (PX) domain, which is a phosphoinositide bi nding domain, and are involved in intracellular trafficking. This protein does not contain a coile d coil region, like some family members, but d oes contain a SH3 domain near its N-terminus. This protein interacts with the cytoplasmic dom ains of the precursor but not the processed for ms of a disintegrin and metalloprotease domain 9 and 15. This protein binds the beta-appendag



	e domain of adaptor protein 2 and may function noto assist adaptor protein 2 in its role at the plasma membrane. This protein interacts with a ctivated Cdc42-associated kinase-2 to regulate the degradation of epidermal growth factor receptor protein.
Applications:	WB
Name of antibody:	SNX9
Immunogen:	Fusion protein of human SNX9
Full name:	sorting nexin 9
Synonyms :	SNX9; MST155; MSTP155; SDP1; SH3PX1; SH3PX D3A; WISP
SwissProt:	Q9Y5X1
WB Predicted band size:	67 kDa
WB Positive control:	HeLa cells
WB Recommended dilution:	500-2000

