

兔抗 SNX11 多克隆抗体

中文名称：兔抗 SNX11 多克隆抗体

英文名称： Anti-SNX11 rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻（-20℃）

宿 主： Rabbit

抗 原： SNX11

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the sorting nexin family. Members of this family contain a phox (PX) domain, which is a phosphoinositide binding domain, and are involved in intracellular trafficking. This protein does not contain a coiled coil region, like some family members. This gene encodes a protein of unknown function. This gene results in two transcript variants differing in the 5' UTR, but encoding the same protein.
Applications:	ELISA, IHC
Name of antibody:	SNX11
Immunogen:	Fusion protein of human SNX11
Full name:	sorting nexin 11

SwissProt:	Q9Y5W9
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and human esophagus cancer
IHC Recommend dilution:	25-100

