

## 兔抗 LBH 多克隆抗体

中文名称: 兔抗 LBH 多克隆抗体

英文名称: Anti-LBH rabbit polyclonal antibody

别名: LBH regulator of WNT signaling pathway

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: LBH

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	Transcriptional activator which may act in mitogen-activated protein kinase signaling pathway.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	LBH
<b>Immunogen:</b>	Fusion protein of human LBH
<b>Full name:</b>	LBH regulator of WNT signaling pathway
<b>SwissProt:</b>	Q53QV2
<b>ELISA Recommended dilutio</b>	5000-10000

<b>n:</b>	
<b>IHC positive control:</b>	Human esophagus cancer
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	12 kDa
<b>WB Positive control:</b>	Rat heart tissue, Mouse skeletal muscle tissue lysates
<b>WB Recommended dilution:</b>	500-2000



