

兔抗 TRIM69 多克隆抗体

中文名称：兔抗 TRIM69 多克隆抗体

英文名称： Anti-TRIM69 rabbit polyclonal antibody

别名： Trif; HSD34; RNF36; HSD-34

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： TRIM69

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the RING-B-box-coiled-coil (RBCC) family and encodes a protein with an N-terminal RING finger motif, a PRY domain and a C-terminal S PRY domain. The mouse ortholog of this gene is specifically expressed in germ cells at the round spermatid stages during spermatogenesis and, when overexpressed, induces apoptosis. Alternatively spliced transcript variants encoding distinct isoforms have been described.
Applications:	ELISA, IHC
Name of antibody:	TRIM69

Immunogen:	Fusion protein of human TRIM69
Full name:	tripartite motif containing 69
Synonyms :	Trif; HSD34; RNF36; HSD-34
SwissProt:	Q86WT6
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human brain
IHC Recommend dilution:	25-100

