

## 兔抗 SNX17 多克隆抗体

中文名称：兔抗 SNX17 多克隆抗体

英文名称：Anti-SNX17 rabbit polyclonal antibody

别名：sorting nexin 17

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SNX17

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the sorting nexin family. Members of this family contain a phosphoinositide binding domain, which is a phosphoinositide binding domain, and are involved in intracellular trafficking. This protein does not contain a coiled coil region, like some family members, but contains a B41 domain. This protein interacts with the cytoplasmic domain of P-selectin, and may function in the intracellular trafficking of P-selectin. Alternative splicing results in multiple transcript variants encoding different isoforms.
<b>Applications:</b>	ELISA, WB, IHC

<b>Name of antibody:</b>	SNX17
<b>Immunogen:</b>	Fusion protein of human SNX17
<b>Full name:</b>	sorting nexin 17
<b>SwissProt:</b>	Q15036
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human colorectal cancer and Human esophagus cancer
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	52 kDa
<b>WB Positive control:</b>	293T cell lysate
<b>WB Recommended dilution:</b>	500-2000



