

兔抗 LCAT 多克隆抗体

中文名称： 兔抗 LCAT 多克隆抗体

英文名称： Anti-LCAT rabbit polyclonal antibody

抗 原： LCAT

储 存： 冷冻（-20℃）

宿 主： Rabbit

反应种属： Human, Mouse, Rat

标 记 物： Unconjugate

相关类别： 一抗

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes the extracellular cholesterol esterifying enzyme, lecithin-cholesterol acyltransferase. The esterification of cholesterol is required for cholesterol transport. Mutations in this gene have been found to cause fish-eye disease as well as LCAT deficiency.
Applications:	ELISA, IHC
Name of antibody:	LCAT
Immunogen:	Fusion protein of human LCAT
Full name:	lecithin-cholesterol acyltransferase
SwissProt:	P04180

ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer
IHC Recommend dilution :	25-100

