

兔抗 ZNHIT6 多克隆抗体

中文名称：兔抗 ZNHIT6 多克隆抗体

英文名称：Anti-ZNHIT6 rabbit polyclonal antibody

别名：BCD1; C1orf181; NY-BR-75

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：ZNHIT6

反应种属：Human Mouse

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	Required for box C/D snoRNAs accumulation involved in snoRNA processing, snoRNA transport to the nucleolus and ribosome biogenesis.
Applications:	WB
Name of antibody:	ZNHIT6
Immunogen:	Synthesized peptide derived from internal of human CA 181.
Full name:	zinc finger, HIT-type containing 6
Synonyms :	BCD1; C1orf181; NY-BR-75
SwissProt:	Q9NWK9
WB Predicted band size:	54 kDa

WB Positive control:	Jurkat cells and Hela cells lysates
WB Recommended dilution:	500-3000

