

兔抗 SNX29 多克隆抗体

中文名称：兔抗 SNX29 多克隆抗体

英文名称： Anti-SNX29 rabbit polyclonal antibody

别名： RUNDC2A; A-388D4.1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SNX29

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

SNX29 also known as RUNDC2A, is a 375 amino protein that contains one RUN domain. RUNDC2 A is encoded by a gene that maps to human chromosome 16p13.13, which encodes over 900 genes and comprises nearly 3% of the human genome. The GAN gene is located on chromosome 16 and, with mutation, may lead to giant axonal neuropathy, a nervous system disorder characterized by increasing malfunction with growth. The rare disorder Rubinstein-Taybi syndrome is also associated with chromosome 16, as is Crohn's disease, which is a gastrointestinal inflammatory condition.

Applications:	ELISA, IHC
Name of antibody:	SNX29
Immunogen:	Synthetic peptide of human SNX29
Full name:	sorting nexin 29
Synonyms :	RUNDC2A; A-388D4.1
SwissProt:	Q8TEQ0
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human liver cancer and Human colon cancer
IHC Recommend dilution:	25-100

