

兔抗 ZNRD2 多克隆抗体

中文名称：兔抗 ZNRD2 多克隆抗体

英文名称： Anti-ZNRD2 rabbit polyclonal antibody

别名： zinc ribbon domain containing 2; p27; SSSCA1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： ZNRD2

反应种属： Human, Mouse


标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This antigen is recognized by a subset of anti-centromere antibodies from patients with scleroderma and/or Sjogren's syndrome. Subcellular localization has not yet been established. [provided by RefSeq, Jul 2008]
Applications:	ELISA, IHC
Name of antibody:	ZNRD2
Immunogen:	Fusion protein of human ZNRD2
Full name:	zinc ribbon domain containing 2
Synonyms:	p27; SSSCA1
SwissProt:	O60232
IHC positive control:	Human tonsil and Human prostate cancer

IHC Recommend dilution:	50-300
-------------------------	--------



mlbio 梅联生物
Good elisakit producers

