

小鼠抗 TIMM29 单克隆抗体

中文名称：小鼠抗 TIMM29 单克隆抗体

英文名称： Anti-TIMM29 mouse monoclonal antibody

别名： translocase of inner mitochondrial membrane 29; TIM29; C19orf52

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： TIMM29

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	Component of the TIM22 complex, a complex that mediates the import and insertion of multi-pass transmembrane proteins into the mitochondrial inner membrane. The TIM22 complex forms a twin-pore translocase that uses the membrane potential as the external driving force. Required for the stability of the TIM22 complex and functions in the assembly of the TIMM22 protein into the TIM22 complex. May facilitate cooperation between TIM22 and TOM complexes by interacting with TOMM40.
Applications:	WB

Name of antibody:	TIMM29
Immunogen:	Fusion protein of human TIMM29
Full name:	translocase of inner mitochondrial membrane 29
Synonyms:	TIM29; C19orf52
SwissProt:	Q9BSF4
WB Predicted band size:	29 KD
WB Positive control:	Human placenta tissue lysate
WB Recommended dilution:	500-2000