

兔抗 ZP1 多克隆抗体

中文名称：兔抗 ZP1 多克隆抗体

英文名称：Anti-ZP1 rabbit polyclonal antibody

别名：OOMD; HEL163

相关类别：一抗

储存：冷冻（-20℃）避光

宿主：Rabbit

抗原：ZP1

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	The zona pellucida is an extracellular matrix that surrounds the oocyte and early embryo. It is composed primarily of three or four glycoproteins with various functions during fertilization and preimplantation development. The protein encoded by this gene ensures the structural integrity of the zona pellucida. Mutations in this gene are a cause of oocyte maturation defect and infertility.
Applications:	WB
Name of antibody:	ZP1
Immunogen:	Synthesized peptide derived from internal of human ZP1.
Full name:	zona pellucida glycoprotein 1 (sperm receptor)

Synonyms :	OOMD; HEL163
SwissProt:	P60852
WB Predicted band size:	70 kDa
WB Positive control:	LOVO cells lysate
WB Recommended dilution:	500-3000

