

兔抗 SOCS5 多克隆抗体

中文名称：兔抗 SOCS5 多克隆抗体

英文名称： Anti-SOCS5 rabbit polyclonal antibody

别名： CIS6; CISH6; Cish5; SOCS-5

抗原： SOCS5

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene contains a SH2 domain and a SOCS BOX domain. The protein thus belongs to the suppressor of cytokine signaling (SOCS) family, also known as STAT-induced STAT inhibitor (SSI) protein family. SOCS family members are known to be cytokine-inducible negative regulators of cytokine signaling. The specific function of this protein has not yet been determined. Two alternatively spliced transcript variants encoding an identical protein have been reported.
Applications:	ELISA, IHC

Name of antibody:	SOCS5
Immunogen:	Fusion protein of human SOCS5
Full name:	suppressor of cytokine signaling 5
Synonyms :	CIS6; CISH6; Cish5; SOCS-5
SwissProt:	O75159
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and human ovarian cancer
IHC Recommend dilution:	25-100

