

小鼠抗 TNFRSF4 单克隆抗体

中文名称：小鼠抗 TNFRSF4 单克隆抗体

英文名称： Anti-TNFRSF4 mouse monoclonal antibody

别名： TNF receptor superfamily member 4; OX40; ACT35; CD134; IMD16; TXGP1L

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Mouse

抗原： TNFRSF4

反应种属： Human

标记物： Unconjugate

克隆类型： mouse monoclonal

技术规格

Background:	The protein encoded by this gene is a member of the TNF-receptor superfamily. This receptor has been shown to activate NF-kappaB through its interaction with adaptor proteins TRAF2 and TRAF5. Knockout studies in mice suggested that this receptor promotes the expression of apoptosis inhibitors BCL2 and BCL2L1/BCL2-XL, and thus suppresses apoptosis. The knockout studies also suggested the roles of this receptor in CD4+ T cell response, as well as in T cell-dependent B cell proliferation and differentiation.
Applications:	WB, FC

Name of antibody:	TNFRSF4
Immunogen:	Fusion protein of human TNFRSF4
Full name:	TNF receptor superfamily member 4
Synonyms:	OX40; ACT35; CD134; IMD16; TXGP1L
SwissProt:	P43489
WB Predicted band size:	29 KD
WB Positive control:	PHA-PBMC cell lysate
WB Recommended dilution:	200-1000
FC Recommend dilution:	50-100