

兔抗 ZPBP 多克隆抗体

中文名称：兔抗 ZPBP 多克隆抗体

英文名称：Anti-ZPBP rabbit polyclonal antibody

别名：ZPBP1

抗原：ZPBP

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	ZPBP is one of several proteins that are thought to participate in secondary binding between acrosome-reacted sperm and the egg-specific extracellular matrix, the zona pellucida (McLeskey et al., 1998 [PubMed 9378618]). May be implicated in gamete interaction during fertilization.
Applications:	ELISA, WB, IHC
Name of antibody:	ZPBP
Immunogen:	Fusion protein of human ZPBP
Full name:	zona pellucida binding protein
Synonyms :	ZPBP1

SwissProt:	Q9BS86
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human brain
IHC Recommend dilution:	25-100
WB Predicted band size:	40 kDa
WB Positive control:	Human testis tissue lysate
WB Recommended dilution:	200-1000



