

## 兔抗 LDHAL6B 多克隆抗体

- 中文名称：兔抗 LDHAL6B 多克隆抗体
- 英文名称：Anti-LDHAL6B rabbit polyclonal antibody
- 别名：lactate dehydrogenase A like 6B; LDHL; LDH6B; LDHAL6
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：LDHAL6B
- 反应种属：Human
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	LDHAL6B belongs to the LDH/MDH superfamily. LDH family. The specific function of the protein remains unknown.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	LDHAL6B
<b>Immunogen:</b>	Fusion protein of human LDHAL6B
<b>Full name:</b>	lactate dehydrogenase A like 6B
<b>Synonyms:</b>	LDHL; LDH6B; LDHAL6
<b>SwissProt:</b>	Q9BYZ2
<b>ELISA Recommended diluti</b>	5000-10000

<b>on:</b>	
<b>IHC positive control:</b>	Human tonsil
<b>IHC Recommend dilution:</b>	50-300

