

兔抗 ZSCAN25 多克隆抗体

中文名称： 兔抗 ZSCAN25 多克隆抗体

英文名称： Anti-ZSCAN25 rabbit polyclonal antibody

别名： ZNF498

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： ZSCAN25

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes a protein that bears some similarity to zinc finger proteins, which are involved in DNA binding and protein-protein interactions. Multiple alternatively spliced transcript variants have been identified, but the full-length nature for most of them has not been determined
Applications:	WB
Name of antibody:	ZSCAN25
Immunogen:	Synthesized peptide derived from internal of human ZNF498.
Full name:	zinc finger and SCAN domain containing 25
Synonyms :	ZNF498

SwissProt:	Q6NSZ9
WB Predicted band size:	61 kDa
WB Positive control:	A549 cells and Huvec cells lysates
WB Recommended dilution:	500-3000

