

## 兔抗 LDHB 多克隆抗体

中文名称： 兔抗 LDHB 多克隆抗体

英文名称： Anti-LDHB rabbit polyclonal antibody

别名： LDH-B; LDH-H; LDHBD; TRG-5

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： LDHB

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

**Background:**

This gene encodes an enzyme which catalyzes the reversible conversion of lactate and pyruvate, and NAD and NADH, in the glycolytic pathway. Mutations in this gene are associated with lactate dehydrogenase B deficiency. Pseudogenes have been identified on the X chromosome and on chromosome. Multiple alternatively spliced variants, encoding the same protein, have been identified. Rare LDHB variants result in deficiency of lactate dehydrogenase, a condition with no deleterious effects on health. LDHB deficiency is of interest to laboratory medicine mainly because it can cause mis

	diagnosis in those disorders in which elevation of serum LDH is expected. Lactate dehydrogenase deficiency can probably be considered a non-disease.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	LDHB
<b>Immunogen:</b>	Synthetic peptide of human LDHB
<b>Full name:</b>	lactate dehydrogenase B
<b>Synonyms :</b>	LDH-B; LDH-H; LDHBD; TRG-5
<b>SwissProt:</b>	P07195
<b>ELISA Recommended dilution:</b>	2000-10000
<b>WB Predicted band size:</b>	37 kDa
<b>WB Positive control:</b>	Hela cells
<b>WB Recommended dilution:</b>	1000-5000

