

## 兔抗 ZSCAN29 多克隆抗体

中文名称：兔抗 ZSCAN29 多克隆抗体

英文名称：Anti-ZSCAN29 rabbit polyclonal antibody

别名：ZNF690; Zfp690

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ZSCAN29

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	May be involved in transcriptional regulation.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	ZSCAN29
<b>Immunogen:</b>	Fusion protein of human ZSCAN29
<b>Full name:</b>	zinc finger and SCAN domain containing 29
<b>Synonyms:</b>	ZNF690; Zfp690
<b>SwissProt:</b>	Q8IWY8
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human prostate cancer and human thyroid cancer
<b>IHC Recommend dilution:</b>	25-100



