

## 兔抗 SOX9 多克隆抗体

中文名称：兔抗 SOX9 多克隆抗体

英文名称：Anti-SOX9 rabbit polyclonal antibody

别名：SRY-box 9; CMD1; SRA1; CMPD1; SRXX2; SRXY10

相关类别：一抗

抗原：SOX9

储存：冷冻（-20℃）

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene recognizes the sequence CCTTGAG along with other members of the HMG-box class DNA-binding proteins. It acts during chondrocyte differentiation and, with steroidogenic factor 1, regulates transcription of the anti-Muellerian hormone (AMH) gene. Deficiencies lead to the skeletal malformation syndrome campomelic dysplasia, frequently with sex reversal.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	SOX9
<b>Immunogen:</b>	Synthetic peptide of human SOX9

<b>Full name:</b>	SRY-box 9
<b>Synonyms:</b>	CMD1; SRA1; CMPD1; SRXX2; SRXY10
<b>SwissProt:</b>	P48436
<b>ELISA Recommended dilution:</b>	5000-10000
<b>WB Predicted band size:</b>	56 kDa
<b>WB Positive control:</b>	293T cell lysate
<b>WB Recommended dilution:</b>	500-2000

