

兔抗 LDLRAD1 多克隆抗体

中文名称：兔抗 LDLRAD1 多克隆抗体

英文名称：Anti-LDLRAD1 rabbit polyclonal antibody

抗 原：LDLRAD1

储 存：冷冻（-20℃）避光

宿 主：Rabbit

相关类别：一抗

反应种属：Human

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	Belongs to the LDLR family. 3 isoforms of the human protein are produced by alternative splicing.
Applications:	WB
Name of antibody:	LDLRAD1
Immunogen:	Synthesized peptide derived from internal of human LDLRAD1.
Full name:	low density lipoprotein receptor class A domain containing 1
Synonyms :	0
SwissProt:	Q5T700
WB Predicted band size:	22 kDa
WB Positive control:	Jurkat cells and HepG2 cells lysates
WB Recommended dilution:	500-3000

