

## 兔抗 LDLRAD2 多克隆抗体

中文名称：兔抗 LDLRAD2 多克隆抗体

英文名称：Anti-LDLRAD2 rabbit polyclonal antibody

抗原：LDLRAD2

储存：冷冻（-20℃）避光

宿主：Rabbit

相关类别：一抗

反应种属：Human Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Belongs to the LDLR family.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	LDLRAD2
<b>Immunogen:</b>	Synthesized peptide derived from internal of human LDLRAD2.
<b>Full name:</b>	low density lipoprotein receptor class A domain containing 2
<b>Synonyms :</b>	0
<b>SwissProt:</b>	Q5SZI1
<b>WB Predicted band size:</b>	29 kDa
<b>WB Positive control:</b>	Jurkat cells, LOVO cells and NIH/3T3 cells lysates
<b>WB Recommended dilution:</b>	500-3000

