

兔抗 LDLRAD3 多克隆抗体

中文名称：兔抗 LDLRAD3 多克隆抗体

英文名称：Anti-LDLRAD3 rabbit polyclonal antibody

别名：LRAD3

抗原：LDLRAD3

储存：冷冻（-20℃）避光

宿主：Rabbit

反应种属：Human Mouse

相关类别：一抗

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	May influence APP processing, resulting in a decrease in sAPP -alpha production and increased amyloidogenic P3 peptide production By similarity.
Applications:	IF
Name of antibody:	LDLRAD3
Immunogen:	Synthesized peptide derived from internal of human LDLRAD3.
Full name:	low density lipoprotein receptor class A domain containing 3
Synonyms :	LRAD3
SwissProt:	Q86YD5
IF positive control:	A549 cells
IF Recommend dilution:	100-500

