

兔抗 ZZZ3 多克隆抗体

中文名称: 兔抗 ZZZ3 多克隆抗体

英文名称: Anti-ZZZ3 rabbit polyclonal antibody

抗 原: ZZZ3

别 名: ATAC1

相关类别: 一抗

储 存: 冷冻(-20℃)

宿 主: Rabbit

反应种属: Human, Mouse

标 记 物: Unconjugate

克隆类型: Unconjugate

技术规格

Background:

ZZZ3 (ZZ-type zinc finger-containing protein 3) is a 903 amino acid protein that contains one HTH myb-type DN A-binding domain and one ZZ-type zinc finger. Phospho rylated upon DNA damage by ATM or ATR, ZZZ3 is a s ubunit of the ATAC complex, which is composed of GC N5, CRP2BP, ADA3, TADA2L, DR1, CCDC101, YEATS2, W DR5 and MBIP. The ATAC complex has histone acetyltra nsferase activity on histones H3 and H4. ZZZ3 is expres sed as four isoforms produced by alternative splicing an d is encoded by a gene mapping to human chromosom e 1. Chromosome 1 is the largest human chromosome s panning about 260 million base pairs and making up 8 % of the human genome. There are about 3,000 genes



	on chromosome 1, and considering the great number of genes there are also a large number of diseases associa ted with chromosome 1. Notably, the rare aging disease Hutchinson-Gilford progeria is associated with the LMNA gene which encodes lamin A. When defective, the LMN A gene product can build up in the nucleus and cause characteristic nuclear blebs. The mechanism of rapidly e nhanced aging is unclear and is a topic of continuing e xploration. Stickler syndrome, Parkinsons, Gaucher diseas e and Usher syndrome are also associated with chromos ome 1.
Applications:	ELISA, IHC
Name of antibody:	ZZZ3
Immunogen:	Fusion protein of human ZZZ3
Full name:	zinc finger, ZZ-type containing 3
Synonyms :	ATAC1
SwissProt:	Q8IYH5
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer
IHC Recommend dilution:	25-100

