

## 兔抗 LEF1 多克隆抗体

中文名称：兔抗 LEF1 多克隆抗体

英文名称：Anti-LEF1 rabbit polyclonal antibody

别名：LEF-1; TCF10; TCF7L3; TCF1ALPHA

抗原：LEF1

储存：冷冻（-20℃）

宿主：Rabbit

相关类别：一抗

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

**Background:**

This gene encodes a transcription factor belonging to a family of proteins that share homology with the high mobility group protein-1. The protein encoded by this gene can bind to a functionally important site in the T-cell receptor-alpha enhancer, thereby conferring maximal enhancer activity. This transcription factor is involved in the Wnt signaling pathway, and it may function in hair cell differentiation and follicle morphogenesis. Mutations in this gene have been found in somatic sebaceous tumors. This gene has also been linked to other cancers, including androgen-independent prostate cancer. Alternative splicing results in multiple transcri

	pt variants.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	LEF1
<b>Immunogen:</b>	Fusion protein of human LEF1
<b>Full name:</b>	lymphoid enhancer-binding factor 1
<b>Synonyms :</b>	LEF-1; TCF10; TCF7L3; TCF1ALPHA
<b>SwissProt:</b>	Q9UJU2
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human cervical cancer and Human tonsil
<b>IHC Recommend dilution:</b>	25-100

