

## 小鼠抗 A4GNT 单克隆抗体

中文名称: 小鼠抗 A4GNT 单克隆抗体

英文名称: Anti-A4GNT mouse monoclonal antibody

别名: alpha4GnT

抗原: A4GNT

储存: 冷冻 (-20℃) 避光

宿主: Mouse

反应种属: Human , Dog , Rat , Monkey , Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein from the glycosyltransferase 3 2 family. The enzyme catalyzes the transfer of N-acetylglucosamine (GlcNAc) to core 2 branched O-glycans. It forms a unique glycan, GlcNAc $\alpha$ 1 $\rightarrow$ 4Gal $\beta$ 1 $\rightarrow$ 6Man $\alpha$ 6 and is largely associated with the Golgi apparatus membrane. [provided by RefSeq, Jul 2008].
<b>Applications:</b>	WB
<b>Name of antibody:</b>	A4GNT
<b>Immunogen:</b>	Fusion protein of human A4GNT
<b>Full name:</b>	alpha-1,4-N-acetylglucosaminyltransferase (A4GNT)
<b>Synonyms:</b>	alpha4GnT

<b>SwissProt:</b>	Q9UNA3
<b>WB Predicted band size:</b>	39 kDa
<b>WB Positive control:</b>	HepG2, Hela, SVT2, A549, COS7, Jurkat, MDCK, PC12, MCF-7 cell lysates
<b>WB Recommended dilution:</b>	500-2000