

## 小鼠抗 ACAA2 单克隆抗体

中文名称: 小鼠抗 ACAA2 单克隆抗体

英文名称: Anti-ACAA2 mouse monoclonal antibody

别名: DSAEC

抗原: ACAA2

储存: 冷冻 (-20°C) 避光

宿主: Mouse

反应种属: Human , Dog , Rat , Monkey , Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: mouse monoclonal

### 技术规格

<b>Background:</b>	The encoded protein catalyzes the last step of the mitochondrial fatty acid beta-oxidation spiral. Unlike most mitochondrial matrix proteins, it contains a non-cleavable amino-terminal targeting signal. [provided by RefSeq, Jul 2008].
<b>Applications:</b>	WB, IHC, IF
<b>Name of antibody:</b>	ACAA2
<b>Immunogen:</b>	Fusion protein of human ACAA2
<b>Full name:</b>	acetyl-CoA acyltransferase 2 (ACAA2)
<b>Synonyms:</b>	DSAEC
<b>SwissProt:</b>	P42765
<b>IHC positive control:</b>	human pancreas tissue and carcinoma of human thyroid tis

	sue; human endometrium tissue and human liver tissue
<b>IHC Recommend dilution:</b>	30-150
<b>WB Predicted band size:</b>	42 kDa
<b>WB Positive control:</b>	HepG2, Hela, SVT2, A549, COS7, Jurkat, MDCK, PC12, MCF-7 cell lysates
<b>WB Recommended dilution:</b>	500-2000