

## 兔抗 LETMD1 多克隆抗体

中文名称: 兔抗 LETMD1 多克隆抗体

英文名称: Anti-LETMD1 rabbit polyclonal antibody

别名: HCCR1; HCCR-1; HCCR-2; HCRR-2; 1110019O13Rik

相关类别: 一抗

储存: 冷冻 (-20℃) 避光

抗原: LETMD1

宿主: Rabbit

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a mitochondrial outer membrane protein. It has a potential role in tumorigenesis, which may result from negative regulation of the p53 tumor suppressor gene. Alternatively spliced transcript variants have been noted for this gene.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	LETMD1
<b>Immunogen:</b>	Fusion protein of human LETMD1
<b>Full name:</b>	LETMD1 domain containing 1
<b>Synonyms :</b>	HCCR1; HCCR-1; HCCR-2; HCRR-2; 1110019O13Rik
<b>SwissProt:</b>	Q6P1Q0

<b>WB Predicted band size:</b>	42 kDa
<b>WB Positive control:</b>	PANC1 cells
<b>WB Recommended dilution:</b>	500-2000

