

## 兔抗 SPAG5 多克隆抗体

中文名称：兔抗 SPAG5 多克隆抗体

英文名称： Anti-SPAG5 rabbit polyclonal antibody

别名： MAP126; DEEPEST; hMAP126

抗原： SPAG5

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

#### Background:

Sperm-associated antigen 5 is a protein that in humans is encoded by the SPAG5 gene. This gene encodes a protein associated with the mitotic spindle apparatus. The encoded protein may be involved in the functional and dynamic regulation of mitotic spindles. Essential component of the mitotic spindle required for normal chromosome segregation and progression into anaphase. Required for chromosome alignment, normal timing of sister chromatid segregation, and maintenance of spindle pole architecture. The astrin (SPAG5)-kinastrin (SKAP) complex promotes stable microtubule-kinetochore attachment.

|                                    |                               |
|------------------------------------|-------------------------------|
|                                    | hments.                       |
| <b>Applications:</b>               | ELISA, IHC                    |
| <b>Name of antibody:</b>           | SPAG5                         |
| <b>Immunogen:</b>                  | Fusion protein of human SPAG5 |
| <b>Full name:</b>                  | sperm associated antigen 5    |
| <b>Synonyms :</b>                  | MAP126; DEEPEST; hMAP126      |
| <b>SwissProt:</b>                  | Q96R06                        |
| <b>ELISA Recommended dilution:</b> | 2000-5000                     |
| <b>IHC positive control:</b>       | Human liver cancer            |
| <b>IHC Recommend dilution:</b>     | 25-100                        |

