

## 兔抗 LGALS4 多克隆抗体

- 中文名称：兔抗 LGALS4 多克隆抗体
- 英文名称：Anti-LGALS4 rabbit polyclonal antibody
- 别名：galectin 4; GAL4; L36LBP
- 相关类别：一抗
- 抗原：LGALS4
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 反应种属：Human
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The galectins are a family of beta-galactoside-binding proteins implicated in modulating cell-cell and cell-matrix interactions. The expression of this gene is restricted to small intestine, colon, and rectum, and it is underexpressed in colorectal cancer.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	LGALS4
<b>Immunogen:</b>	Synthetic peptide of human LGALS4
<b>Full name:</b>	galectin 4
<b>Synonyms:</b>	GAL4; L36LBP
<b>SwissProt:</b>	P56470

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human ovarian cancer and Human thyroid cancer
<b>IHC Recommend dilution:</b>	30-150
<b>WB Predicted band size:</b>	36 kDa
<b>WB Positive control:</b>	HT-29 cell lysate
<b>WB Recommended dilution:</b>	500-2000



