

兔抗 SPARCL1 多克隆抗体

中文名称：兔抗 SPARCL1 多克隆抗体

英文名称：Anti-SPARCL1 rabbit polyclonal antibody

别名：SC1; MAST9; PIG33; MAST 9

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SPARCL1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	SPARCL1 also known as high endothelial venule protein (Hevin) or MAST9, is a 664 amino acid member of the SPARC family of proteins. Highly expressed in lymph node, heart, lung, brain, skeletal muscle, ovary, colon and small intestine, SPARCL1 is a secreted protein that contains one EF-hand domain, one follistatin-like
Applications:	ELISA, IHC
Name of antibody:	SPARCL1
Immunogen:	Synthetic peptide of human SPARCL1
Full name:	SPARC-like 1 (hevin)

Synonyms :	SC1; MAST9; PIG33; MAST 9
SwissProt:	Q14515
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and Human breast cancer
IHC Recommend dilution:	25-100

