

兔抗 LGALS1 多克隆抗体

中文名称：兔抗 LGALS1 多克隆抗体

英文名称：Anti-LGALS1 rabbit polyclonal antibody

别名：galectin like; GRP; HSPC159

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：LGALS1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Does not bind lactose, and may not bind carbohydrates.
Applications:	ELISA, IHC
Name of antibody:	LGALS1
Immunogen:	Synthetic peptide of human LGALS1
Full name:	galectin like
Synonyms:	GRP; HSPC159
SwissProt:	Q3ZCW2
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human ovarian cancer and Human thyroid cancer
IHC Recommend dilution:	50-300

