

小鼠抗 ACVRL1 单克隆抗体

中文名称: 小鼠抗 ACVRL1 单克隆抗体

英文名称: Anti-ACVRL1 mouse monoclonal antibody

别名: activin A receptor like type 1; HHT; ALK1; HHT2; ORW2; SKR3; ALK-1; TSR-I; ACVRLK1

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Mouse

抗原: ACVRL1

反应种属: Human

标记物: Unconjugate

克隆类型: mouse monoclonal

技术规格

Background:	This gene encodes a type I cell-surface receptor for the TGF-beta superfamily of ligands. It shares with other type I receptors a high degree of similarity in serine-threonine kinase subdomains, a glycine- and serine-rich region (called the GS domain) preceding the kinase domain, and a short C-terminal tail. The encoded protein, sometimes termed ALK1, shares similar domain structures with other closely related ALK or activin receptor-like kinase proteins that form a subfamily of receptor serine/threonine kinases. Mutations in this gene are associated with hemorrhagic telangiectasia type 2, also known as Rendu-Osler-Weber syndrome 2.
--------------------	--

Applications:	WB
Name of antibody:	ACVRL1
Immunogen:	Fusion protein of human ACVRL1
Full name:	activin A receptor like type 1
Synonyms:	HHT; ALK1; HHT2; ORW2; SKR3; ALK-1; TSR-I; ACVRLK1
SwissProt:	P37023
WB Predicted band size:	56 KD
WB Positive control:	Jurkat, MCF7, K562 and Hela cell lysates
WB Recommended dilution:	500-2000